

SD GlucoNavi™ Blood Glucose Testing Meter

With the SD GlucoNavi™ Blood Glucose Monitoring System you will find quick results with a small sample size. GlucoNavi will provide pre and post meal markers and Hypo warnings.

- No Coding - Reduces errors from improper coding.
- Gold Coated Test Strips - Accurate, precise results
- AST (Alternate Site Testing) - Capture blood from other sites other than the finger tip
- 5 Second Results - Provides fast, accurate results in only 5 seconds
- Small Blood Sample - Makes it less painful to the patient.
- Results Range 20-600 - Expanded glucose range for reliable results
- Glucose Oxidase Chemistry - Eliminates cross reactivities and interferences
- Large Memory - Store 400 results with date and time
- Hypoglycemic Warning - User defined level (40, 50, 60 or 70 mg/dl)
- Alarm Clock - 4 user programmable alarms
- Average - 7, 14, 30 - day pre/post meal or total averages


- Test method - Glucose Oxidase biosensor
- Blood amount required - 0.9ul
- Blood reference type - Plasma reference
- Measuring unit - mg/dl, mmol/L
- Battery - 3V battery (CR2032)
- SD (Diabetes Management Software) for easy diabetes data management.
- High tech wide pure gold electrode and laser patterning



“GOLD” Advantage

Gold Electrodes -

Wide gold electrodes for more accurate, precise results



Laser Engraved Patterning - Ensures excellent precision

Fast Draw Technology

- Less error due to insufficient sample

A. MGM1	SD GlucoNavi Meter	\$9.00
B. MGS50	SD GlucoNavi Test Strips, 50/pkg.	\$10.00
MGCS1	SD GlucoNavi Control Solution (High and Normal)	\$8.00
C. ML25G200	Lancet, 25 ga., 200 pkg.	\$26.00
D. ML28G200	Lancet, 28 ga., 200/pkg.	\$26.00
E. ML30G200	Lancet, 30 ga., 200/pkg.	\$26.00

Call: 866-231-1222

www.maxpertmedical.com



Expert Products for Expert Care

