




**Expert Products For
Expert Care®**



**Simple, creative and trusted medication management and
medication safety product solutions.**



Designed with a dynamic crushing and grinding action to easily deliver finely crushed medication from the top down.

- Easy dual action crushing and grinding.
- Exclusive top down crushing into bottom cup. Cup provides convenience for mixing and drawing into a feeding syringe.
- Crush directly into a 30mL plastic med cup.
- Free spin, pepper grinder action, delivers fine powder.
- Safe! Latex-free and FDA compliant materials.
- Dimensions: 3-3/8"L x 2-3/16"W
- Sold individually.
- White Blue: BMG02WB
- White Red: BMG02WR

Latex Free

*Patented

THE FIRST OF ITS KIND hand-held twist pill crusher that allows medication to be delivered in two ways:

Ridges on the grinding surface provide fast, dual action crushing and grinding with minimal residue remaining.

Contact us for large quantity pricing.

1-19 20-59 60-119 120+

BMG02WB White Blue \$7.99 \$6.55 \$5.20 \$4.25

BMG02WR White Red \$7.99 \$6.55 \$5.20 \$4.25

Directly into a 30 mL disposable plastic medication cup



Into the bottom snap on clear cup




Comfortable and easy to grip symmetrical design – Ergonomic design and easy to grip top provide comfort when crushing and aid in the twisting effort.



Free Spin – Once the top is screwed down to initially break and crush pills, a simple reverse twist then allows the grinding surface to freely spin similar to a pepper grinder, continuing to grind medications without any resistance and with total ease.

THE FIRST OF ITS KIND hand-held twist pill crusher that allows medication to be delivered in two ways:



Two Delivery Options– Fine powder can be delivered directly into the bottom serving cup where the caregiver empties powder directly into food or adds water and draws into a syringe. The second option is inserting a 30mL medication cup in the crusher well and allowing medication to be crushed directly into the med cup. This minimizes transfer and cleanup.

\$9.99 each
Contact us for quantity and system pricing!




- Dimensions: 3-3/8"L x 2-3/16"W
- Sold individually.
- BMG01B - Blue, BMG01R - Red

Comfortable and easy to grip anti-slip symmetrical design – Ergonomic design and easy to grip rubber material provide comfort when crushing and aid in the twisting effort.



*Patented



Safer, more complete delivery of medication!
Ridges on the grinding surface provide fast, dual action crushing and grinding with minimal residue remaining. Less work for the caregiver and more assurance of medication delivery.



Top Pill Storage Compartment with Hinged Lid - Great for storing pills when traveling.

Internal Cleaning Feature
If you continue to free spin after crushing is complete, air circulates within the crusher blowing remaining powder into the bottom cup, resulting in minimal medication remains on the grinding surface and in the crushing well.



- Easy dual action crushing and grinding for the caregiver delivers a fine powder with minimal residue remaining. Excellent for mixing and feeding tube use.
- Comfortable anti-slip rubber grip.
- Exclusive top down crushing into serving cup.
- Free spin, pepper grinder action, delivers fine powder.
- Safe! Latex-free and FDA compliant materials.
- Built-in storage with hinged lid.

What do you call an expertly designed pill splitter that delivers on quality and function... and includes an exclusive built-in pill puncher feature? We call it the MAXSPLIT Pill Splitter with Puncher proudly crafted by Maxpert Medical. This patent pending splitter is one of a kind.

Quality and durable FDA complaint food grade plastic ensures that medications are contacting high quality, safe plastic and that the splitter is going to last over repeated use. When opening MAXSPLIT, you will find an ergonomic U-shaped holder that allows pills of various sizes, shapes and thicknesses to be easily positioned and split safely and accurately. This U-shaped opening assists in proper pill alignment and the rubber insert holds pill in place while splitting to help ensure an accurate cut.

Have you ever seen a splitter with a pill puncher tool built-in? You have now. This feature can only be found in the MAXSPLIT Pill Splitter with Puncher.

Clear Blue



Clear Red

		1-19	20 - 59	60 - 119	120 +
BMS01B	Clear Blue	\$6.99	\$5.75	\$4.55	\$3.75
BMS01R	White Red	\$6.99	\$5.75	\$4.55	\$3.75

Contact us for large quantity pricing.

Pill Puncher

Inside the splitter's lid is a precision, sharp stainless-steel blade designed to cut through with minimal chances of crumbling or splintering medication. Simple to use! The outside of the lid features a non-slip press point lined area to help identify where to apply light pressure when splitting. Only light pressure is required due to the quality sharp blade. Portable, and compact! MAXSPLIT is scrub or coat pocket friendly and includes a convenient storage compartment for split pills or other medications that require splitting.



- Safe, high quality materials! Latex free, FDA compliant food grade material that is very durable.
- Easy pill placement! The U-shaped holder allows pills of various sizes, shapes and thicknesses to be easily positioned and split safely and accurately. The area is lined with a rubber insert to help pill stay put during splitting.
- Only light pressure required! Stainless-steel precision blade easily cuts through pills with reduced chance of splintering and crumbling.

- Transparent plastic provides an easy view while splitting. Choose your splitter with a Blue or Red rubber insert and compartment cover. BMS01B or BMS01R
- **Bonus pill puncher feature!** This keyhole is an integrated pill puncher feature that is a simple, clever and convenient way to assist in pill removal from blister packs or unit dose packaging.
- Easy to clean with warm water. Rinse and air dry.



RUBBER
STAY PUT
INSERT



Safety is not only delivered with the integrated blade guard to prevent accidental blade contact, but also with the premium FDA complaint food grade materials that MAXSPLIT is produced; latex free, transparent plastic. You can feel safe that MAXSPLIT will protect you from blade contact and confident to know that your medications are in contact only with high quality, safe plastic.

When opening MAXSPLIT, you will find an ergonomic U-shaped holder that allows pills of various sizes, shapes and thicknesses to be easily positioned and split safely and accurately. This U-shaped opening with alignment line assists in proper pill alignment and the rubber insert holds the pill in place while splitting to help ensure an accurate cut.

Clear Red



Clear Blue

- User safety! The integrated blade guard conceals the sharp blade, preventing accidental contact and providing a safe pill splitting experience.
- Safe, high quality materials! Latex free, FDA compliant food grade material that is very durable.
- Easy pill placement! The U-shaped holder allows pills of various sizes, shapes and thicknesses to be easily positioned and split safely and accurately. The area is lined with a rubber insert that helps the pill to stay put while splitting.
- Only light pressure required! Stainless-steel precision blade easily cuts through pills with reduced chance of splintering and crumbling.
- Transparent plastic provides an easy view while splitting. Choose your splitter with a Blue or Red rubber insert and compartment cover. BMS02B or BMS02R

Have you ever seen a splitter with a pill puncher tool built-in? You have now. This feature can only be found in the MAXSPLIT Safety Pill Splitter with Puncher.

- **Bonus pill puncher feature!** This keyhole is an integrated pill puncher feature that is a simple, clever and convenient way to assist in pill removal from blister packs or unit dose packaging.
- Easy to clean with warm water. Rinse and air dry.




PILL PUNCHER



		Each
BMS02B	Clear Blue	\$8.99
BMS02R	White Red	\$8.99

Contact us for large quantity pricing.



Blade Guard
Conceals
Blade

Inside the splitter's lid, safely and discretely hidden by the blade guard is a precision, sharp stainless- steel blade designed to cut through with minimal chances of crumbling or splintering medication. The unique blade guard travels with the blade and keeps it contained to prevent accidental blade contact.

Simple to use! The outside of the lid features a non-slip press point lined area to help identify where to apply light pressure when splitting. Only light pressure is required due to the quality sharp blade.

Portable, and compact! MAXSPLIT is travel friendly and includes a convenient storage compartment for split pills and travel medications that require splitting. Its durable construction will hold up nicely during travel and everyday use.

FOLD THE FLAPS FOR CRUSHING SAFETY

GREAT FOR USP <800> - CRUSHING HAZARDOUS MEDICATIONS

MAXCRUSH® Pill Crusher Pouch

- This specially designed, disposable pill crusher pouch effectively eliminates the problem of cross-contamination. The opening of the pouch forms two flaps. One folds in to provide a pouch barrier, containing contents and the other is used for easy gripping while crushing and then can be folded over after crushing to create a natural closure, containing crushed medication.
- The pouch shape permits easy opening and dispensing of crushed medication and its overall thickness and durability withstand repetitive crushing.
- Universal to pouch pill crushers.
- Clear Polyethylene; Latex free
- 3.9" x 2.0"
- Patent pending
- Sold 1,000 per box; 50 pouches per bag and 20 bags per case



Fold the flap to provide a barrier that minimizes the risk of powder residue being emitted during crushing.

MAXCRUSH®

A. BCB2803 MAXCRUSH Pill Crusher Pouch, 1,000 \$19.99

MAXCRUSH™ M200 Cup Pill Crusher

This expert designed cup pill crusher combines an ergonomic handle design with a unique buffering and grinding motion to deliver highly effective pill crushing with minimal manual effort.

- A safe, simple and effective pill crushing solution.
- Very easy to use. The unique anvil design working with the cup platform delivers both a crushing and grinding with maximum crushing effort. No tapping or pounding.
- Delivers a very fine powder that is ideal for consumption in food or drink or for feeding tube use.
- A two cup method is used to reduce the chance of cross-contamination and particulate powder. Accommodates a range of cup sizes. The 1 oz. (30mL) plastic cups are recommended.
- Delivers a silent and powerful crushing experience.
- The buffering design provides user assistance and absorbs crushing force, thereby reducing chances of wrist and hand related stress injuries.
- MAXCRUSH Cup Pill Crusher is portable, lightweight and has a compact space saving design which is convenient for medication carts, nursing stations, in-room use and near automated dispensing cabinets.
- This crusher is used in various areas of healthcare when the end result desired is quiet, pain free operation.
- Specifications: 1.8 lbs.
- 9-3/4"L x 3-7/8"W x 5-1/2"H



The unique anvil design working with the cup platform delivers both a crushing and grinding with maximum crushing effort. No tapping or pounding.

B. BMC2802 MAXCRUSH Cup Pill Crusher \$149.00

Supply Bin Indicator Dial

The Supply Bin Indicator Dial is a simple and cost-effective way to communicate the supply level status of bins, shelves, cabinets, refrigerators and more. By providing a simple visual alert, a stock out can be prevented and the appropriate inventory level can be maintained.

This Supply Bin Indicator Dial can help prevent wasted money when supplies aren't available, and the time personnel spends trying to find out whether supplies have been re-ordered.

Simply peel and stick! Remove the film from the adhesive on the backside and attach to your desired area.

A simple turn of the indicator dial from Green to Red and Red to Yellow provides visual communication of supply status. Green = Supply level is good, Red = Supply level is low and it's time to reorder, Yellow = Supplies have been ordered.

- Self-adhesive or magnetic mount provides versatility and mounting convenience.
- Simple to use and a great complement to existing systems such as automation.
- Colors in the Supply Bin Indicator Dial have simple meaning and help effectively communicate status.
- The affordability of this visual supply management tool allows it to be introduced across multiple areas within a department or supply chain. The Supply Bin Indicator Dial can effectively communicate supply status in medication rooms, the central pharmacy, materials management, supply rooms, warehouses and many more areas.
- Blue plastic with labeling area. The large area can be written on with a marker or labeled, making it convenient to number bins and boxes for further management.
- Available individually.


SpinFlag and BinFlag Supply Level Indicators

The SpinFlag and BinFlag help save time, money and effort with your supply management system. They help with First in First Out processes and monitoring of stock levels, as well as alerts for restocking. These indicators help provide a simple alert before stock has run out, helping to maintain the appropriate stock levels.

Attached to open bins, shelves, cabinets etc. as a simple method to identify when important supplies are in need of restocking or have been ordered for replenishment.

- Extremely simple installation as the back plate supplied can be fixed to any type of bin, cabinet, shelf or fixture.
- Very visual, creating a reduction in stock checking time.
- Easy to use and compliments your current supply management system.
- Can be affixed with glue, screws, double sided tape or other methods. Each indicator has a fixture plate that mounts and then the indicator slides into the plate.
- Manufactured in hi-impact resistant ABS plastic making them durable and robust, perfect for any tough environment.
- SpinFlag Dimensions: 1.4" x 1.4"
Small BinFlag Dimensions: 1.7"L x 3.5"W
Large BinFlag Dimensions: 1.9"L x 4.4"W
- Sold individually

A. DS104	Supply Bin Indicator Dial	\$2.95 each
B. DS101	SpinFlag Stock Indicator	\$8.95 each
C. DS102	BinFlag Stock Indicator, Small	\$10.95 each
D. DS103	BinFlag Stock Indicator, Large	\$11.50 each




SpinFlag



BinFlag



Adjustable Pull Tight Security Seals




- This adjustable pull tight security seal or drug lock provides tamper indication and security for crash carts, crash cart trays, med kits, medical bags and boxes while providing a quick visual indication when something has been opened or used.
- Consecutive number option assists with documentation and control.
- This one-piece seal is designed to be easily applied and removed. It requires no tools for removal, saving valuable time when accessing medications or supplies.
- This all plastic seal has a unique engineered break point that allows quick access in emergencies.


A.	DTE03B	Adjustable Pull Tight Security Seal, Numbered, Blue, 100/pkg.	\$12.00
B.	DTE03G	Adjustable Pull Tight Security Seal, Numbered, Green, 100/pkg.	\$12.00
C.	DTE03R	Adjustable Pull Tight Security Seal, Numbered, Red, 100/pkg.	\$12.00
D.	DTE03Y	Adjustable Pull Tight Security Seal, Numbered, Yellow, 100/pkg.	\$12.00
E.	DTE04B	Adjustable Pull Tight Security Seal, No Numbers, Blue, 100/pkg.	\$12.00
F.	DTE04G	Adjustable Pull Tight Security Seal, No Numbers, Green, 100/pkg.	\$12.00
G.	DTE04R	Adjustable Pull Tight Security Seal, No Numbers, Red, 100/pkg.	\$12.00
H.	DTE04Y	Adjustable Pull Tight Security Seal, No Numbers, Yellow, 100/pkg.	\$12.00

- Thin seal tail fits various applications such as zipper ring, latches, eyelets and holes.
- Adjustable Pull Tight Security Seals are flexible by design; engineered to bend and not break when applying.
- Average breaking strength: 10-15 lbs.
- Seal length: 9"
- Usable length: 6-1/2"
- Polypropylene plastic
- 10 seals per strip; 100 seals per package
- Available in Blue, Green, Red and Yellow.

Expiration Pull Tight Security Seal




- Large expiration print with white write on block allows an expiration date to be written in with a pen or marker; thereby eliminating the need to open and search through items to identify expiration dates or apply additional expiration indicators such as labels.
- Consecutive numbers for documented control
- Can easily be applied and removed. Requires no tools for removal, saving valuable time accessing important supplies.
- Thin adjustable seal tail fits numerous rings, zippers, latches and holes.
- Average breaking strength - 11-14 lbs.
- Seal length - 8" - Tail length - 7"
- Autoclavable up to 325° F.
- Available in Blue, Green, Red and Yellow.



- Designed to indicate whether an item has been opened or compromised.
- Features low breaking strength; break without cutting.
- Simple, one-piece polypropylene 9"L seal. Provides tamper-evident protection for boxes, cases, cabinets, fire exit doors, fire extinguishers, first aid boxes and more.
- For maximum security, these barcoded seals offer a unique sequential numbering system that cannot be rubbed off.

- Average breaking strength - 10-15 lbs.
- Code 128 barcode print
- Available in Orange.
- Packed 100 per package.

A.	DTE01B	Expiration Seal, Blue, 100/pkg.	\$12.00
B.	DTE01G	Expiration Seal, Green, 100/pkg.	\$12.00
C.	DTE01R	Expiration Seal, Red, 100/pkg.	\$12.00
D.	DTE01Y	Expiration Seal, Yellow, 100/pkg.	\$12.00
E.	DTE02	Adjustable Barcode Seal, Orange	\$12.00



Adjustable Barcode Pull Tight Security Seal, Numbered

Steri-Tamp® Sterile Tamper Evident Seals


- The only single-use tamper-evident seal that provides a 100% sterile barrier.
- 3X stronger adhesion than other seals on the market
- Protects against contamination and accidental double dosing.
- Alerts staff when a drug has been added.
- Patented dual-layer leaves true tamper evidence with “OPENED” labeling left behind.
- The foil top layer provides evidence of needle puncture.
- Allows pharmacy to recycle medication and protects patients from dosing errors.
- Maintains adhesion, even when refrigerated or frozen in cold drug storage conditions (down to -20 degrees centigrade).
- Packaged in easy-to-use dispensers, containing 1,000 units on a roll.



Steri-Tamp® Vial Seals

- Notches allow vial seals to lay flat to prevent bacteria from entering the port.
- Also protects Pharmacy Technicians from Cytotoxins

Item #	Size in mm	Color	Quantity	Price
DTE09	13 mm	Red	1000/ roll	\$125.00
DTE08	20 mm	Silver	1000/ roll	\$125.00
DTE07	28 mm	Blue	1000/ roll	\$125.00



Steri-Tamp® IV Bag Port Seals

- Protects the point of entry of IV admixture bags from contamination and accidental double dosing.
- Using a simple “Twist” method, the seal wraps around the port and itself.

Item #	Color	Quantity	Price
DTE05	Blue	1000/ roll	\$125.00
DTE06	Red	1000/ roll	\$125.00
DTE12	Yellow	1000/ roll	\$125.00



Steri-Tamp® Belly Button Port Seal

- Covers the belly button port specifically of a Hospira IV bag and can also be used on 13mm vials.
- Stronger adhesion made to work on the Belly Button Bag & vials.

Item #	Size in mm	Color	Quantity	Price
DTE11	13 mm	Green	1000/ roll	\$125.00



Steri-Tamp® Non-Sterile Tamper Evident Seals

Steri-Tamp® Red Syringe Seal

- Non-sterile solution to provide tamper evidence for syringes and other medical containers.
- Alerts staff when the syringe has been opened or used.
- Improves workflow; saves time and money.
- Helps prevent risk of double dosing.
- Allows pharmacy to recycle medication.
- Maintains adhesion, even when refrigerated or frozen in cold drug storage conditions (down to -20 degrees centigrade).
- Packaged in easy-to-use dispensers, containing 1000 units each.
- Smart choice for USP compliance.

Item #	Color	Quantity	Price
DTE10	Red	1000/ roll	\$125.00



Steri-Tamp® Tamper-Clear Syringe Seal™

- Clear visibility to the barrel of the syringe.
- Useful for oral syringes, epi-pens, inhalers and other devices with markings that need to be seen or scanned.
- Non-sterile solution to provide tamper evidence for syringes and other medical containers.
- Improves workflow; saves time and money.
- Alerts staff when the syringe has been opened or used.
- Helps prevent risk of double dosing.
- Allows pharmacy to recycle medication.
- Maintains adhesion, even when refrigerated or frozen in cold drug storage conditions (down to -20 degrees centigrade)
- Packaged in easy-to-use dispensers, containing 1000 units per roll. Smart choice for USP compliance

Item #	Color	Quantity	Price
DTE13	Red Tint	1000/ roll	\$170.00



Tamper Resistant Tape

- Excellent deterrent to tampering!
- This affordable tape provides immediate visual detection of tampering.
- After application tamper resistant tape is designed to tear easily into multiple pieces, even if removed immediately. The serrated edges (saw tooth pattern) enhance the security by making the tape break off into irregular little strips.
- The words Tamper Resistant Tape appear in a continuous pattern on the tape roll.
- Translucent red acetate material with black print.
- Ideal for bottles, jars, vials, tubes, envelopes, boxes, bags and more.
- Each cardboard dispensing box includes two 0.5” x 216 ft (72 yds) rolls.

A. DTE16	Tamper Resistant Tape, Red 0.5 in x 216 ft., box of 2 rolls	\$38.00
----------	--	---------

Security Tamper Tape

- Yellow is ideal for sealing hazardous or chemotherapy bags, vials, cases, and more as a deterrent to tampering!
- Security Tamper Tape features write on sections for Sealed By, Date and has a white write on strip which provide additional customizing of security features. Can be written on with permanent marker or ink pen.
- This affordable tape provides immediate visual detection of tampering.
- After application security tamper tape is designed to tear easily into multiple pieces, even if removed immediately. The serrated edges (saw tooth pattern) enhance the security by making the tape break off into irregular little strips.
- The words Security Seal repeats every five inches throughout the roll.
- This tape can be used as a continuous tape or cut or torn into custom lengths.
- Translucent yellow acetate material with black print.
- The paper backing has a crack-and-peel feature for convenience when removing from the roll.
- Ideal for bottles, jars, vials, tubes, envelopes, boxes, bags and more.
- Each cardboard dispensing box includes one 1.25” x 108 ft roll.

B. DTE17	Security Tamper Tape, Yellow 1.25 in. x 108 ft., 1 roll	\$42.00
----------	--	---------

TAMPER RESISTANT TAPE



Crash Cart/Drug Tray Security Bags

- These tamper evident crash cart tray bags help ensure the integrity and security of medications and supplies.
- Bags are a convenient and easy means to provide tamper evidence to crash cart and anesthesia trays. No lids or seals to mess with or manage.
- Bags are simple to seal by simply removing the liner across the strong adhesive strip and folding the lip of the bag over.
- Crash Cart/Drug Tray Security Bags allow trays to be prepared ahead of time and available for quick exchange when a crash cart is opened or case is completed.
- An integrated perforation in the center of the bag is highlighted for quick identification and provides quick and easy access to the tray contents in an emergency.
- Available in 2 common sizes that accommodate the majority of crash cart tray sizes. 14” x 22” ideal for half size style trays or boxes and 20” x 29” geared toward the full size trays or boxes.
- 2 mil thick
- Available in convenient package sizes of 20 bags per package.




A. DTE14	Crash Cart/Drug Tray Security Bag, 14 X 22, 20/pkg.	\$5.00
B. DTE15	Crash Cart/Drug Tray Security Bag, 20 X 29, 20/pkg.	\$10.50



Tamper Evident Transport Bags

- These tamper evident bags help maintain integrity of medications and supplies during transport.
- Once bags are sealed they must be torn to open, providing clear indication of use. The tape closure is aggressive, causing the top of the bag to tear and wrinkle up if there is an attempt to open from the top.
- Tamper Evident Bags feature a center perforation for quick, convenient access when opening. 2 mil thick. Available in 6”W x 6”H, 8”W x 10”H and 10”W x 14”H
- Tamper Evident Tear Receipt Bags feature a white write on block and perforated tear receipt with write on block so the bag can be written on with identical information to safeguard against tampering. These bags are torn at the top perforation to open after being sealed. 1.5 mil thick. Available in 3.75”W x 5.75”H and 6”W x 10”H.
- Available 100 per package.

C. DTE18	Tamper Evident Bag, 6”W x 6”H, 2mil, 100/pkg.	\$18.00
D. DTE19	Tamper Evident Bag, 8”W x 10”H, 2mil, 100/pkg.	\$24.00
E. DTE20	Tamper Evident Bag, 10”W x 14”H, 2mil, 100/pkg.	\$32.00
F. DTE21	Tamper Evident Tear Receipt Bag, 3.75”W x 5.75”H, 1.5mil, 100/pkg.	\$7.95
G. DTE22	Tamper Evident Tear Receipt Bag, 6”W x 10”H, 1.5mil, 100/pkg.	\$12.00



Maxpert Tamper Evident Reusable Bags



- Provide security and maintain integrity of your medications, supplies, and other important items during handling and transport with these durable and sealable tamper evident reusable bags.
- Maxpert Tamper Evident Reusable Bags are used in conjunction with pull tight or padlock style security seals to make them tamper evident. Eyelets on the high quality two zipper system interlock allowing a seal to be inserted, preventing the zippers and bag from being opened without breaking the seal.
- Two styles are available. All Clear No Mesh and Clear with Mesh
- Clear with mesh are high quality, clear bags made of a mesh-embedded transparent material which allows the user to identify if the material has been cut, indicating tampering has occurred. The embedded mesh also enhances the bag strength for security and longevity of the bag.
- Enable tracking and control by utilizing a security seal with consecutive numbering. Use with Maxpert Medical's DTE03 Numbered Pull Tight Security Seal or any of our other pull tight seals. No tools required to remove our breakable pull tight seals.

A. DTE29	Tamper Evident Reusable Bag, No Mesh, 11" x 7"	\$8.50
B. DTE30	Tamper Evident Reusable Bag, With Mesh, 11" x 7"	\$8.50
C. DTE31	Tamper Evident Reusable Bag, No Mesh, 13" x 9"	\$9.50
D. DTE32	Tamper Evident Reusable Bag, With Mesh, 13" x 9"	\$9.50

High Alert Medication Transport Bags

These reusable printed bags with high alert medication message are designed to safely and securely transport high alert or look-alike sound-alike medications.

- These bags segregate high alert medications needing transported from the pharmacy to medication rooms, automated dispensing cabinets, patient rooms and alternate care sites.
- Red backing allow these bags to stand out and the medical professional to take notice that the medications enclosed are high alert and should be handled appropriately with double checks.
- Grommet in upper corner allows a security seal to be inserted through the grommet and zipper pull to secure the bag. A broken seal may indicate tampering or use.
- Reusable! Durable vinyl material can easily be wiped clean and withstand multiple uses.
- Narrow Bag: 10-1/2"W x 3-1/2"H, Large Bag: 11"W x 9"H

E. AMS02	Narrow Bag	\$7.50
F. AMS03	Large Bag	\$8.50



- Reusable means you help eliminate waste and save per use versus tamper evident disposable bags.
- Bag fronts feature an external access window for important content identification, department id or routing information. Window is 4"W x 3-1/8"H
- Expanded bottom permits easier loading and adds a little extra storage space for transported items.
- Constructed of durable PVC material.



Patient MDI Bag

- Whether the metered dose inhaler is dispensed in the hospital or the patient brings their own from home, the MDI needs a convenient way to be identified and kept clean. This pre-printed bag with bold lettering alerts that the bag contains an MDI and is great for transporting from the pharmacy to the patient's room, to SPD for sterilizing, or back to pharmacy for cleaning and potential reuse.
- The bag is convenient for keeping track of an MDI in the patient room whether it's their own or being used during admission.
- Bag is clear with white block and blue print.
- Ideal for transporting and clean storage.
- Available in 6" X 9", 2 mil, 100 per package



Do Not Return To Pharmacy Bag

- Utilize this pre-printed reclosable bag to transport, store and collect medication that is going into or has been in a patient isolation room or an isolation area. The bold message not to return medication to the pharmacy or automated dispensing cabinet serves as a cautionary message to ensure overall safety. This simple warning bag helps ensure that if medications have been exposed to contamination, they do not get transported back to the pharmacy or entered bag into an automated dispensing cabinet where cross contamination can quickly spread.
- Bag is clear with white block and bold orange print.
- Great for isolation precautions such as COVID-19!
- Ideal for transporting syringes, vials, IV bags and related supplies.
- Available in 6" X 9", 2 mil, 100 per package

A. EBG02B	Patient MDI Bag, 100/pkg.	\$19.90
B. EBG01B	Do Not Return To Pharmacy Bag, 100/pkg.	\$19.90

Opioid Overdose Kit Bag

- Use this durable laminated nylon opioid overdose kit bag for assembling your Naloxone kits! Ample storage for kit supplies and large white print on the exterior, Opioid Overdose Kit, clearly identifies the bag and contents to be used for prevention of opioid drug overdoses.
- Dimensions: 11" x 6"



C. AMS01	Opioid Overdose Kit Bag	\$7.95
D. AMS37	Hazardous Drug Spill Kit Supply Bag	\$32.00

Hazardous Drug Spill Kit Supply Bag

- This durable, bright yellow supply bag allows you to keep multiple hazardous drug spill kits and PPE organized, easily identified and readily available to respond to hazardous drug breakages and spills.
- Dimensions: 18" x 14" x 4"





Maxpert NIOSH Listed HD Transport Bag

Improve pharmacy operation procedures and your USP 800 safety initiatives with these durable and sealable bags dedicated to safe NIOSH Listed Hazardous Drug transport.

- These heavy-duty thick bags feature a dedicated hazardous drug caution message along with bold color and warning to alert users to safely handle and observe cautions that go along with handling and transport of hazardous medications.
- Bags are resealable and leak resistant.
- Ideal for transporting syringes, vials, IV bags and related supplies from the pharmacy to med or patient rooms.
- Value priced for long term savings.
- Available in four sizes: 4" x 6" 4mil, 6" x 9" 4mil, 9" x 12" 4mil, 12" x 15" 4mil
- Sold 100 per package.

Maxpert Clear Reclosable Bag w/ Yellow Block

- Available in four sizes for color coding: 3" x 5", 4" x 6", 6" x 9", 9" x 12"
- 4mil thick. Sold 100 per package.



AMSH246	NIOSH Listed HD Transport Bag, 4"W x 6"H, 4mil, 100/pkg.	\$9.99
AMSH269	NIOSH Listed HD Transport Bag, 6"W x 9"H, 4mil, 100/pkg.	\$10.99
AMSH2912	NIOSH Listed HD Transport Bag, 9"W x 12"H, 4mil, 100/pkg.	\$19.99
AMSH21215	NIOSH Listed HD Transport Bag, 12"W x 15"H, 4mil, 100/pkg.	\$29.99
AMS29	Clear Reclosable Bag w/ Yellow Block, 3" x 5", 4mil, 100/pkg.	\$9.50
AMS24	Clear Reclosable Bag w/ Yellow Block, 4" x 6", 4mil, 100/pkg.	\$9.99
AMS25	Clear Reclosable Bag w/ Yellow Block, 6" x 9", 4mil, 100/pkg.	\$10.99
AMS26	Clear Reclosable Bag w/ Yellow Block, 9" x 12", 4mil, 100/pkg.	\$19.99

Maxpert Caution Hazardous Drug Transport Bags

- These heavy-duty thick bags feature a dedicated caution hazardous drug message along with bold color and warning to alert users to safely handle and observe cautions that go along with handling, disposal and transport of hazardous medications.
- Bags are resealable and leak resistant.
- Ideal for transporting syringes, vials, IV bags and related supplies from the pharmacy to med or patient rooms.
- Value priced for long term savings.
- Available in four sizes: 4" x 6" 4mil, 6" x 9" 4mil, 9" x 12" 4mil, 12" x 15" 4mil
- Sold 100 per package.

Improve pharmacy operation procedures and your USP 800 safety initiatives with these durable and sealable bags dedicated to Hazardous Drug transport.



AMSH146	Caution Hazardous Drug Transport Bag, 4"W x 6"H, 4mil, 100/pkg.	\$9.99
AMSH169	Caution Hazardous Drug Transport Bag, 6"W x 9"H, 4mil, 100/pkg.	\$10.99
AMSH1912	Caution Hazardous Drug Transport Bag, 9"W x 12"H, 4mil, 100/pkg.	\$19.99
AMSH11215	Caution Hazardous Drug Transport Bag, 12"W x 15"H, 4mil, 100/pkg.	\$29.99

Help comply with Section II of USP 800



A



B



C



D



E

Chemotherapy Drug Transport Bags

- These extra heavy-duty thick bags feature a dedicated chemotherapy caution message along with color and symbol to alert users to safely handle and observe cautions that go along with chemotherapy and hazardous medications.
- Bags are resealable and leak resistant.
- Ideal for transporting cytotoxic medications and related supplies from the pharmacy to med or patient rooms.
- 4 mil thick plastic
- Value priced for long term savings.
- Available in four sizes: 6" x 9", 9" x 12", 12" x 15", 13" x 19"
- Sold 100 per package.

A.	AMS04	Chemo Drug Transport Bag, 6"W x 9"H, 4mil, 100/pkg.	\$10.99
B.	AMS05	Chemo Drug Transport Bag, 9"W x 12"H, 4mil, 100/pkg.	\$19.99
C.	AMS06	Chemo Drug Transport Bag, 12"W x 15"H, 4mil, 100/pkg.	\$29.99
D.	AMS07	Chemo Drug Transport Bag, 13"W x 19"H, 4mil, 100/pkg.	\$43.99
E.	AMS08	Chemo Drug Transport Bag (yellow), 12"W x 15"H, 4mil, 100/pkg.	\$29.99

USP 800 Hazardous Drug and Spill Kit Labels

NIOSH Listed HD:
Follow facility handling precautions.

HAZARDOUS DRUGS

USP 800
HAZARDOUS DRUG

USP 800
HAZARDOUS DRUG
OBSERVE SPECIAL HANDLING, ADMINISTRATION & DISPOSAL REQUIREMENTS
1 x 3

USP 800
HAZARDOUS DRUG
OBSERVE SPECIAL HANDLING, ADMINISTRATION & DISPOSAL REQUIREMENTS
1 x 4

AMS28	NIOSH Listed HD Label, 1-9/16" x 5/16", 1,000/bx	\$6.00
AMS33	Red Hazardous Drugs Label, 1-9/16" x 5/16", 1,000/bx	\$6.00
AMS22	USP 800 HD Label 1-9/16" x 5/16", 1,000 per box	\$6.00
AMS20	USP 800 HD Label 3" x 1", 1,000 per box	\$12.00
AMS21	USP 800 HD Label 4" x 1", 1,000 per box	\$14.00
AMS32	Hazardous Drug Spill Kit Label, 6" x 4-1/4, each	\$1.75



Hazardous Drug Safety Caution Signs

CAUTION

OBSERVE SAFETY PRECAUTIONS

HAZARDOUS DRUGS ARE HANDLED AND STORED IN THIS AREA

CAUTION

OBSERVE SAFETY PRECAUTIONS

HAZARDOUS DRUGS ARE HANDLED AND STORED IN THIS AREA

HD

CAUTION

OBSERVE SAFETY PRECAUTIONS

HAZARDOUS DRUGS ARE HANDLED AND STORED IN THIS AREA

HD

SPILL KIT

HD

SPILL KIT

- Signs are printed on durable 3mm PVC and include easy removable strips for mounting.
- Contact us for your custom sign needs.

Hazardous Drug Caution Window Clings

CAUTION

OBSERVE SAFETY PRECAUTIONS

HAZARDOUS DRUGS ARE HANDLED AND STORED IN THIS AREA

HD

CAUTION

OBSERVE SAFETY PRECAUTIONS

HAZARDOUS DRUGS ARE HANDLED AND STORED IN THIS AREA

- Durable window clings apply like a label and can easily be removed and repositioned. (not reverse cling)

Hazardous Drug Caution Magnets

CAUTION

HAZARDOUS DRUG

OBSERVE SAFETY PRECAUTIONS FOR HANDLING AND ADMINISTRATION

CAUTION

HAZARDOUS DRUG

OBSERVE SAFETY PRECAUTIONS FOR HANDLING AND ADMINISTRATION

CAUTION

HAZARDOUS DRUG

OBSERVE SAFETY PRECAUTIONS FOR HANDLING AND ADMINISTRATION

AMS14	Caution Hazardous Drug Sign, 10" x 7"	\$6.00
AMS27	Caution Hazardous Drug Sign, 6" x 4-1/4"	\$2.00
AMS16	Hazardous Drug Sign with Triangle, 10" x 7"	\$6.00
AMS30	Hazardous Drug Spill Kit Sign, 10" x 7"	\$6.00
AMS31	Hazardous Drug Spill Kit Sign, 6" x 4-1/4"	\$2.00
AMS13	Caution Hazardous Drug Window Cling, 10" x 7"	\$8.00
AMS15	Hazardous Drug Window Cling with Triangle, 10" x 7"	\$8.00
AMS18	Caution Hazardous Drug Magnet, 2.5" x 1.5"	\$0.70
AMS19	Caution Hazardous Drug Magnet, 5" x 1.5"	\$0.75
AMS19	Caution Hazardous Drug Magnet, 5" x 1.5"	\$0.75

PPE Reminder Sign, Label and Window Cling

Safety When Using Automated Dispensing Cabinets

PPE REMINDER FOR AUTOMATED DISPENSING CABINET USERS

- REMOVE GLOVES
- PERFORM HAND HYGIENE
- ACCESS MACHINE & REMOVE MEDS
- WIPE TOUCH POINTS WITH DESIGNATED WIPES
- PERFORM HAND HYGIENE
- DON NEW GLOVES IF APPROPRIATE

PPE REMINDER FOR AUTOMATED DISPENSING CABINET USERS

- REMOVE GLOVES
- PERFORM HAND HYGIENE
- ACCESS MACHINE & REMOVE MEDS
- WIPE TOUCH POINTS WITH DESIGNATED WIPES
- PERFORM HAND HYGIENE
- DON NEW GLOVES IF APPROPRIATE

Surface contamination is a huge concern with COVID-19 and in general. There is a risk of transferring the virus to the ADC keyboard as well as any meds that are touched and for surfaces beyond the automated dispensing cabinet.

Simply post a sign, label or window cling as a reminder for proper protocol and PPE use when using automated dispensing cabinets.

AMS34	PPE Reminder Sign, 6" x 4-1/4"	\$2.00
AMS35	PPE Reminder Label, Removable Vinyl, 6" x 4-1/4"	\$1.75
AMS36	PPE Reminder Window Cling, 6" x 4-1/4"	\$2.00

Contains Tranexamic Acid Warning Labels

ISMP recently released a NAN alert regarding wrong-route errors with tranexamic acid. One of the listed error prevention measures is to place a Contains Tranexamic Acid label on the vial.

Contains Tranexamic Acid

Contains Tranexamic Acid

AMS38	White/Blk	\$6.00 pkg.
AMS39	Red/White	\$6.00 pkg.

- 1,000 labels per package
- 1-1/2"W x 1/2"H



DRUG DISPOSAL
FAST • EASY • SAFE
FOR THE ENVIRONMENT



Pharmacist Approved & Recommended

Rx Destroyer™ Instant Drug Disposal System

HOW Rx Destroyer™ WORKS...

The Rx Destroyer™ is eco-friendly “disposal in a bottle”. The Rx Destroyer™ bottle contains a ready-to-use chemical digestion solution. Simply load unused or unwanted medications into the bottle. Shake to mix patent pending digestion solution over the medications. Drugs start dissolving on contact permitting the active medication ingredients to be adsorbed and neutralized by the active charcoal ingredients. The active medication ingredients permanently bond to the activated charcoal and the process is irreversible.

Rx Destroyer™ is reusable and is recommended to be stored in a controlled area when not in use. Discard bottle with contents into common trash.

REGULATIONS

Rx Destroyer meets “non-retrievable” DEA standards for controlled substances. In addition, Rx Destroyer meets DEA requirements for drug destruction which helps accomplish the DEA's goal to reduce drug diversion by keeping medications out of the wrong hands, diminish accidental poisoning and decrease illicit drug activities. **Available in All-Purpose or Liquids Only.**

Rx Destroyer™ is the responsible drug disposal choice and it supports EPA and FDA efforts to help prevent leaching of pharmaceutical into drinking water, soil, and waterways. In most states, Rx Destroyer™ meets disposal regulations for solid waste, which allows the Rx Destroyer™ bottle and its contents to be discarded with common trash. For additional topic related information, please refer to your facilities, state and local regulations.



- Patent Pending Ventilated Cap
- Safe to use and safe for the environment
- Accepts all non-hazardous medications including narcotics
- Makes medications irreversible



The 64 oz. and 1 Gallon Pro Series All-Purpose and Liquids Only Version has a tethered cap to prevent loss and view gauge to monitor the fill level so you know when full.



Available in a universal all purpose or liquids only version.



- 5 gallon for liquids holds approx 500 oz., 2.5 gallon for liquids holds 250 oz. of liquids and the 64 oz. gallon for liquids holds 50 oz. of liquids.

AMD01	4 oz. bottle (holds 50 pills)	\$4.40
AMD03	16 oz. bottle (holds 300 pills)	\$9.25
AMD02	64 oz. bottle, 4/cs. (holds 1500 pills)	\$87.00
AMD18	Pro-Series All-Purpose 64 oz. bottle, 4/cs.	\$99.00
AMD04	64 oz. bottle, liquids only, 4/cs. (includes 4 packs of hardener)	\$121.00
AMD19	Pro-Series Liquids Only 64 oz. bottle, 4/cs. (includes hardener packets)	\$133.00
AMD07	1 gallon Pro All-Purpose, 4/cs. (holds 3000 pills)	\$205.00
AMD08	1 gallon Pro Liquids Only, 4/cs. (holds 100 oz. of liquids)	\$233.00
AMD05	2.5 gallon bottle, 2/cs. (holds 7500 pills)	\$209.00
AMD06	2.5 gallon bottle, liquids only, 2/cs. (includes 2 packs of hardener)	\$241.00
AMD12	5 gallon All Purpose bottle (holds 15,000 pills)	\$216.00
AMD13	5 gallon Liquids only bottle (holds 500 oz.)	\$247.00
AMD09	Rx Destroyer Bottle Funnel	\$7.00
AMD10	Wall Mount Bracket for 64 oz. Bottle	\$48.00
AMD11	Wall Mount Bracket for 64 oz. Bottle	\$19.00